

The versatile Cyto system from Hettich is a valuable support for applications in cytological diagnostics and offers a number of advantages:

- **Maximum safety: Cyto rotor with certified biocontainment**

Cytologists have been waiting for it for a long time, and now it's finally here: the cyto rotor that can be sealed against biological hazards. The New Hettich cyto rotor 1515-A is the first to have certified biocontainment*. This means maximum safety and optimum protection against infectious Samples in centrifuging.

*) Tested by the TÜV in conformity with DIN EN 61010, section 2-020

- **Time savings:** With the new 12-place rotors up to 12 samples can be centrifuged in a single run.

- **SUBSTANTIAL COST SAVINGS**

Cytology laboratories that have used alternative products to date (such as Shandon® centrifuges) can now use just a single Hettich centrifuge. There is no need to acquire a centrifuge solely for cytology purposes. Users who are already Hettich customers can use the new rotors in their existing ROTOFIX 32 A and UNIVERSAL 320 models.

- **MAXIMUM FLEXIBILITY**

The 12-place cyto are easy to insert and remove without the need for any tools.

