

## SCHOOL LABORATORY TABLE



Available in 1000mm, 1200mm and 1500mm widths. The height is 900mm; choice of levelling feet or heavy duty braked castors.

A fully welded 38mm x 38mm x 1.5mm box section construction with a central bottom rail. All metalwork is cleaned and etched being sprayed with a powder coat paint (colour options: White/Grey/Black) that is oven baked to an attractive and durable finish. This design is capable of heavy duty static loading.

Fitted with a 16mm Trespa TopLab worktop (colour options: White/Speckled White/Grey/ Speckled Grey/Silver) All edges are chamfered and polished. Trespa TopLab is chemically robust providing excellent scratch and abrasion resistance. It is suitable for temperatures ranging from -40°C to + 180°C.