



Heating Mantles

Thermal Insulation:	Moulded silica powder
Heating Element:	Nichrome wire
Mantle:	Woven glass fibre mat
Element	Temperature: Up to 450°C (Max)
Power Supply:	220V AC

Stirred	Volume (ml)	Load (W)	Unstirred
	50	60	E-1005N
ES-301N	100	80	E-101N
ES-302N	250	150	E-102N
ES-303N	500	260	E-103N
ES-304N	1000	350	E-104N
ES-305N	2000	450	E-105N
ES-306N	3000	550	E-106N
ES-307N	5000	700	E-107N
ES-308N	10000	1200	E-108N
ES309N	20000	1800	E-109N
O140UK	Support rod for any of the above		