

Clifton Unstirred Water Baths

NE1D Series Unstirred Digital Water Baths

Clifton Range®



NE1D - 14



Comprising of a high quality corrosion resistant stainless steel tank, with a bright clean finish housed in a durable painted outer case. Painted surfaces feature anti-bacterial finish with good chemical resistance. The NE1D series, Clifton Digital Baths are designed for immersion of flasks/bottles/racks in a highly stable temperature environment throughout the water bath.

	NE1D-4	NE1D-8	NE1D-14	NE1D-22
GENERAL				
Temp. Range	Ambient + 5°C - 99°C			
Temp. Control	Digital, LED display showing actual or set temperature, resolution 1°C			
Sensitivity/Uniformity	±0.2°C / ±0.1°C			
Heater Power (220-240V)	400W	800W	1000W	1500W
Voltage	120V or 230V models available 50-60Hz			
Drain	N/A			Drain Outlet
DIMENSIONS				
Tank Capacity	4 Litre	8 Litre	14 Litre	22 Litre
External (mm)	332w x 185d x 290h	332w x 270d x 290h	361w x 332d x 290h	538w x 332d x 290h
Internal (mm)	300w x 150d x 150h	300w x 240d x 150h	325w x 300d x 150h	500w x 300d x 150h
Water Level Max. (mm)	130			
Water Level Min. (mm)	55			

NE1D Series Unstirred Digital Water Baths

Accessories - Lids



Stainless Steel Gabled Lid with Insulated Handle

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Accessory Code	SL1-4	SL1-8	SL1-14	SL1-22



Stainless Steel Flat One Piece Lid with Insulated Handle

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Accessory Code	LD-4	LD-8	LD-14	LD-22



Stainless Steel Hinged Gabled Lid with Insulated Handle

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Accessory Code	N/A		SL1-14H	SL1-22H



Clear Hinged Gable Lid

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Accessory Code	GL1-4	GL1-8	GL1-14	GL1-22



Clear Lift Off Gable Lid

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Accessory Code	GL1-4AF	GL1-8AF	GL1-14AF	GL1-22AF



Stainless Steel Lid with Concentric Rings

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Accessory Code	SLR1-4	SLR1-8 SLR2-8	SLR1-14 SLR2-14 SLR3-14	SLR1-22 SLR3-22

SLR1-4: 2 x 114mm holes

SLR1-8: 4 x 83mm holes

SLR2-8: 2 x 114mm holes

SLR1-14: 4 x 114mm holes

SLR2-14: 6 x 83mm holes

SLR3-14: 4 x 114mm holes

SLR1-22: 6 x 114mm holes

SLR3-22: 6 x 114mm holes



Polypropylene Spheres (Pack of 200)

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
No. of packs	1 Pack	2 Packs		4 Packs
Accessory Code	BP0368			


NE1D Series Unstirred Digital Water Baths

Accessories



Stainless Steel Test Tube Racks

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
No. of Racks	1	2	4	6
Accessory Code	6870 : 26 holes x 17mm DIA			
	6871 : 16 holes x 26mm DIA			
	6872 : 36 holes x 13mm DIA			
	6873 : 18 holes x 19mm DIA			
	6900 : 12 holes x 32mm DIA			
	6920 : 9 holes x 40mm DIA			



Stainless Steel Flat n Fold™ Test Tube Racks

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
No. of Racks	1	2	4	6
Accessory Code	6770 : 26 holes x 17mm DIA			
	6771 : 16 holes x 26mm DIA			
	6772 : 36 holes x 13mm DIA			
	6773 : 18 holes x 19mm DIA			
	6800 : 14 holes x 32mm DIA			



Churchill Test Tube Racks

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Racks Available	C18-4 C25-4	C18-8 C25-8	C18-14 C25-14	N/A
Accessory Code	C18-4 : 28 holes x 18mm DIA			
	C18-8 : 44 holes x 18mm DIA			
	C18-14 : 68 holes x 18mm DIA			
	C25-4 : 28 holes x 25mm DIA			
	C25-8 : 44 holes x 25mm DIA			
	C25-14 : 68 holes x 25mm DIA			



Drain Tap

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Accessory Code	N/A			BX0953

NE1D Series Unstirred Digital Water Baths

Accessories



Stainless Steel Wire Basket

Bath	NE1D-4	NE1D-8	NE1D-14	NE1D-22
Accessory Code	BSK-4	BSK-8	BSK-14	BSK-22

Universal Accessories

TC-1	Thermometer Clip complete with Bent Stem Spirit Filled Thermometer
CAL	Calibration Record
BX0616	Draining Syphon
BKF-1	Bar Keepers Friend for cleaning water baths and removing dirt
LB-2.5	2.5 Litres Lab Bath 4590 Heat Transfer Fluid
LB-5.0	5.0 Litres Lab Bath 4590 Heat Transfer Fluid