

# Clifton Stirred Water Baths

## NE4 - D/CT Series Stirred Digital Baths

*Clifton Range®*



NE4-22D/CT

Clifton NE4D/CT digital stirred baths come complete with sturdy clear polycarbonate tank, allowing observation of samples. Complete with stainless steel shelf designed for good water circulation and temperature uniformity. All Clifton stirred baths feature a low liquid level float switch, which ensures heating is switched off if the liquid falls below the recommended level.

### NE4 - 14D/CT

### NE4 - 22D/CT

GENERAL	
Temp. Range	Ambient + 5°C - 90°C (-0°C with accessory cooling)
Temp. Alarm	Under/over temperature alarm is automatically set 4°C above and 4°C below set point
Temp. Control	PID Digital, LED display, resolution 0.1°C
Sensitivity/Uniformity	±0.05°C / ±0.01°C
Timer	Visual indicator and end message - Variable 00:01 to 99:59 hh:mm
Heater Power (220-240V)	1250W
Voltage	120V or 230V models available 50-60Hz
DIMENSIONS	
Tank Capacity	14 Litre
External (mm)	530w x 325d x 295h
Internal (mm)	385w x 290d x 150h
Water Level Max	130

# NE4-D/CT Series Stirred Digital Baths

## Accessories



### Polypropylene Spheres (Pack of 200)

Bath	NE4-14D/CT	NE4-22D/CT
No. of Packs	2 Packs	3 Packs
Accessory Code	BP0368	



### Immersion Dip Cooler [ambient to -20°C] \*

Bath	NE4-14D/CT	NE4-22D/CT
Accessory Code	DC1-300	

\* Please note : NE4D/CT range will only operate to 0°C with accessory cooling.




### Cooling Coil

Bath	NE4-14D/CT	NE4-22D/CT
Accessory Code	NE4-CC	



### Stainless Steel Test Tube Racks

Bath	NE4-14D/CT	NE4-22D/CT
No. of Racks	3	5
Accessory Code	6870 : 26 holes x 17mm DIA	
	6871 : 16 holes x 26mm DIA	
	6872 : 36 holes x 13mm DIA	
	6873 : 18 holes x 19mm DIA	
	6900 : 12 holes x 32mm DIA	
	6920 : 9 holes x 40mm DIA	



### Stainless Steel Flat n Fold™ Test Tube Racks

Bath	NE4-14D/CT	NE4-22D/CT
No. of Racks	3	5
Accessory Code	6770 : 26 holes x 17mm DIA	
	6771 : 16 holes x 26mm DIA	
	6772 : 36 holes x 13mm DIA	
	6773 : 18 holes x 19mm DIA	
	6800 : 14 holes x 32mm DIA	

## Universal Accessories

CAL	Calibration Record
BX0616	Draining Syphon
LB-2.5	2.5 Litres Lab Bath 4590 Heat Transfer Fluid
LB-5.0	5.0 Litres Lab Bath 4590 Heat Transfer Fluid