

E SERIES –40°C UPRIGHT FREEZERS

MODEL	V080-40E	V145-40E	V255-40E	V495-40E
EXT DIMS (mm) H x W x D	720 x 555 x 670	940 x 600 x 670	1540 x 600 x 735	1940 x 750 x 790
CAPACITY (litres)	80	145	255	495
DRAWERS	2	4	6	14*

- Adjustable temperature range –25°C to –40°C
- Electronic control
- Digital temperature display
- Audio/visual alarm system
- Temperature deviation and mains fail alert
- Back up temperature indicator
- Volt free alarm contacts
- Solid pull out drawers
- Extra thick insulation
- Manual defrost
- Options: Chart recorder, castors, door lock, access port for monitoring probe
- V495-40E fitted with lock and castors as standard

V080-40E



V495-40E



* V495-40E can be fitted with 3 inner sub doors instead of drawers