

K SERIES ECONOMY -40°C UPRIGHT FREEZERS

MODEL	V080-40K	V145-40K	V255-40K	V495-40K
EXT DIMS (mm) H x W x D	630 x 555 x 670	850 x 600 x 670	1450 x 600 x 735	1850 x 750 x 790
CAPACITY (litres)	80	145	255	495
DRAWERS	2	4	6	14*

- Adjustable temperature range -25°C to -40°C
- Internally sparkfree
- Solid pull out drawers
- Manual defrost
- Electronic control
- Digital temperature display
- Visual temperature alarm
- Extra thick insulation

Options: Audible temperature alarm, access port for monitoring probe, chart recorder, castors

* V495-40K can be fitted with 3 inner sub doors instead of drawers — specify when ordering

V495-40K



V255-40K



V080-40K

