



REFRIGERATION TECHNOLOGY

RLE RANGE OF +4°C REFRIGERATORS

MODEL	SW200 RLE	SW400 RLE	SW550 RLE	SW600 RLE
EXT DIMS (mm) H x W x D	850 x 600 x 600	1850 x 600 x 600	1720 x 780 x 720	1890 x 780 x 735
CAPACITY (litres)	200	400	550	600
SHELVES	3	4	4	4
WEIGHT (kgs)	45	75	85	90

- Adjustable temperature range +2°C to +8°C
- Electronic control
- Digital temperature display
- Fan circulation for even temperature distribution
- Visual temperature alarm
- Adjustable wire shelves
- Lockable door
- Automatic defrost
- White exterior
- Optional stainless steel exterior
- Glass door available on SW200 RLE
- Options: Audible alarm, access port for monitoring probe
- These models are not internally sparkfree

SW400 RLE



SW200 RLE (with glass door)



PETER SWAN & SONS LTD, 3 DRYDEN LOAN, BILSTON GLEN IND ESTATE, LOANHEAD EH20 9HR

TEL: 0131 448 0880

www.peterswan.co.uk

FAX: 0131 448 0881