

ULT –85°C UPRIGHT FREEZERS

MODEL	ULT100-U	ULT415-U	ULT580-U	ULT725-U
EXT DIMS (mm) H x W x D	830 x 900 x 530	1930 x 800 x 855	1940 x 1025 x 855	1950 x 1025 x 830
CAPACITY (litres)	100	415	580	725
COMPARTMENTS	2	4	5	5
WEIGHT (kgs)	115	235	265	315

- Adjustable temperature range –55°C to –85°C
- Electronic control
- Digital temperature display
- Audio/visual alarm system
- Temperature deviation and mains fail alert
- Fault diagnosis facility
- Volt free alarm contacts
- Adjustable shelves (fixed on ULT100-U)
- Lockable door
- Inner sub doors
- Stainless steel interior
- Manual defrost
- Access port for monitoring probe
- Options: Chart recorder, CO2 back-up, Storage racks

ULT580-U



ULT100-U

