

ULT –85°C CHEST FREEZERS

MODEL	ULT580-C	ULT650-C	ULT760-C
EXT DIMS (mm) H x W x D	1090 x 1690 x 790	1075 x 2080 x 850	1090 x 2050 x 790
CAPACITY (litres)	580	650	760
WEIGHT (kgs)	240	280	395

- Adjustable temperature range –55°C to –85°C
- Electronic control
- Digital temperature display
- Audio/visual alarm system
- Temperature deviation and mains fail alert
- Fault diagnosis facility
- Volt free alarm contacts
- Lockable lid
- Inner sub lids
- Stainless steel interior
- Manual defrost
- Access port for monitoring probe
- Options: Chart recorder, CO2 back-up, Storage racks

Models ULT580-C and ULT760-C have vacuum insulation panels

ULT580-C



ULT760-C



PETER SWAN & SONS LTD, 3 DRYDEN LOAN, BILSTON GLEN IND ESTATE, LOANHEAD EH20 9HR

TEL: 0131 448 0880

www.peterswan.co.uk

FAX: 0131 448 0881