

VFD RANGE OF -25°C UPRIGHT FREEZERS

MODEL	SW80 VFD	SW130 VFD	SW220 VFD
EXT DIMS (mm) H x W x D	630 x 555 x 625	850 x 555 x 625	1450 x 600 x 630
CAPACITY (litres)	80	130	220
DRAWERS	2	3	6

MODEL	SW260 VFD	SW360 VFD	SW520 VFD
EXT DIMS (mm) H x W x D	1650 x 600 x 630	1750 x 700 x 750	1725 x 750 x 750
CAPACITY (litres)	260	360	520
DRAWERS	7	7	14

- Adjustable temperature range -15°C to -25°C
- Internally sparkfree
- Solid pull out drawers
- Adjustable thermostat
- Temperature alarm
- Manual defrost

Options: Electronic temperature control, audible alarm, door lock, castors, access port for monitoring probe



SW80 VFD

SW130 VFD

SW220 VFD

SW520 VFD